

# Quick Set Up Guide For S-6000

Thank you for your service!

This **Sound Oasis®** unit was designed to help you manage tinnitus and improve your sleep.

It includes 288 sound options, an AM/FM/clock/radio/calendar and Bluetooth® technology so you can stream music/sounds/APPS from your Bluetooth device. A dedicated S-6000 APP provides even more tinnitus therapy options.


## Quick Set Up Steps

- 1 Press the SOUNDS button to play sounds. Press [▲] or [▼] until you reach your favorite sound. Rotate the VOLUME/ON/OFF dial to adjust sound volume.
- 2 Set the time by pressing the SET TIME button for 2 seconds. Press [▲] or [▼] to reach your desired time/date and the SET TIME button to select.
- 3 Set the alarm time/type by pressing the SET 1 button for 2 seconds. Press [▲] or [▼] to reach your desired alarm time/sound and the SET 1 button again to select.
- 4 Press STEREO button to play AM/FM radio. Select AM/FM with the AM/FM button. Press [▲] or [▼] until you reach your favorite radio station.

- 5 To turn on the Bluetooth feature, press the SOUNDS button (*the last played sound will play*). You then have 5 minutes to pair your Bluetooth enabled device with your S-6000. Select S-6000 from the device list on your Bluetooth enabled device (*the time and BLUETOOTH will be shown on the S-6000 display when pairing is successful*). Once you start playing sound/music from your Bluetooth enabled device, it will start playing from your S-6000.
- 

## The S-6000 APP includes:

- 48 'made for tinnitus' sounds
- a 12-band equalizer so you can customize sound frequency levels.



Go to Apple's APP Store  
or Google Play  
to download



Sound Oasis is passionate about helping you manage your tinnitus. If you have any questions or product service issues, please contact us at: **1-866-625-3218** or email us at **info@soundoasis.com** | **www.soundoasis.com**