Massager
Massagegeräte
Masseur
Masajeador
Massaggiatore
Masażer
Massageapparaat

JS10036

USA office: Fontana **AUS** office: Truganina **ITA** office: Milano **POL** office: Gdańsk
DEU office: FDS GmbH, Neuer Höltingbaum 36, 22143 Hamburg, Deutschland **FRA** office: Saint Vigor d'Ymonville
GBR office: FDS Corporation Limited, Unit 4, Blackacre Road, Great Blakenham, Ipswich, Ip6 0FL, United Kingdom

THIS INSTRUCTION BOOKLET CONTAINS **IMPORTANT** SAFETY INFORMATION.
PLEASE READ AND KEEP FOR FUTURE REFERENCE.

EN

Before You Start

- ⚠ Please read all instructions carefully.
- ⚠ Retain instructions for future reference.
- ⚠ Separate and count all parts and hardware.
- ⚠ Read through each step carefully and follow the proper order.
- ⚠ We recommend that, where possible, all items are assembled near to the area in which they will be placed in use, to avoid moving the product unnecessarily once assembled.
- ⚠ Always place the product on a flat, steady and stable surface.
- ⚠ Keep all small parts and packaging materials for this product away from babies and children as they potentially pose a serious choking hazard.

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How to Use:

EN

1. Fully charge the massager before first use. Only use the charger provided and do not use the device while charging. The power button flashes while charging and will have a solid green light when fully charged.
2. Choose the preferred massage head from the 3 heads included in the package. Simply align the required head shaft with the socket and firmly press until the shaft is fully inserted. To remove, simply pull the attachment to withdraw from the socket.
3. Press the power button for 3 seconds to turn on the massager. The default mode is massage mode 1 and speed intensity level 2 to adjust the settings, simply press the massage mode and intensity level buttons.
4. The massager is pre-programmed to run for 15 minutes and will automatically shut off after the session. If additional session is desired, simply restart the massager.

Note: The massager has 3 interchangeable massage heads that can be used in different parts of the body:

- A. Flat massage head for shaping of various muscle parts.
- B. Double ball head for massaging the acupoints on the body.
- C. W-shaped massage head for massaging arms, waist, back buttocks, thighs.

How to Determine the Battery Capacity of the Massager:

When the unit is off, short press the power button and check the light indicator:

Flashing Red Light: 0-5%, (unit won't turn on)

Steady Red Light for 5 seconds: 6-20% (a prompt to charge the unit)

Steady Yellow Light for 5 seconds: 21-50%

Steady Green Light for 5 seconds: More than 50%

When the unit is on, just look at the power button and the color display will show the battery capacity.

Red Light: Under 20%

Yellow Light: 20-50%

Green Light: Over 50%

Product Overview

• Flat Massage Head

• Double Ball Head

• PU Self-Formed
Leather Head

• Charging Cord

Power Button: Press the power button for 3 seconds to start the massager.

Massage Mode Button: When you turn on the massager, it automatically starts in Massage Mode 1. You can switch between the 3 modes by pressing the Massage Mode Button.

Intensity Level Button: The default intensity level is Level 2. Press the Intensity Level Button to switch the intensity levels.

Note: The massager has 3 interchangeable massage heads that can be used in different parts of the body:

- A. Flat massage head for shaping of various muscle parts.
- B. Double ball head for massaging the acupoints on the body.
- C. W-shaped massage head for massaging arms, waist, back buttocks, and thighs.

Troubleshooting

PROBLEM	POSSIBLE CAUSE	SOLUTION
My unit suddenly stopped working	Power outage OR massage more than 15min OR long working time, the unit has been shut off due to overheating protection	The unit automatically shuts off after 15min of use. Use the unit after the motor has been naturally cooled
Unit stopped working and will not turn on after pressing the power button.	The unit is out of power or the battery shuts off due to overcurrent protection.	If the unit is out of power, charge it again before using. If it still has battery and will not turn on, charge it for a few seconds to activate.
Unit stopped charging	Unit not receiving the correct amount of voltage	Make sure you use the charger included in the package. If issue persists, contact us
For more troubleshooting guides, please scan the QR code.		

WARNING! CAUTION! DANGER! FIRE HAZARD!

CHARGING OR DISCHARGING BATTERIES POSES A NUMBER OF SIGNIFICANT HAZARDS WHICH CAN INCREASE THE RISK OF FIRE, EXPLOSION, SERIOUS INJURY, AND PROPERTY DAMAGE IF NOT PERFORMED AS PER THE INSTRUCTIONS.

BEFORE USE, READ AND FOLLOW ALL MANUFACTURER'S INSTRUCTIONS, WARNINGS AND PRECAUTIONS.

NEVER ALLOW CHILDREN UNDER THE AGE OF 14 TO USE THE MASSAGER WITHOUT THE SUPERVISION OF A RESPONSIBLE, KNOWLEDGEABLE ADULT.

THIS IS NOT A TOY.

Care, Storage & Cleaning:

After use, wipe the massager with a soft, damp cloth using non-abrasive cleaners and store it in a cool and dry place.

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The following batteries or accumulators are contained in this electronic device

EN

Battery Type	Chemical System
3S1P2000	Li-ion

Information on the safe removal of the batteries or accumulators

- Warning: Make sure the battery is completely drained.
- Remove the battery in the power-off state.
- Carefully remove the battery or accumulator.
- The battery or accumulator and the device can now be disposed of separately.

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EN

Return / Damage Claim Instructions

- ⚠ **DO NOT discard the box / original packaging.**
In case a return is required, the item must be returned in original box. Without this your return will not be accepted.
- ⚠ **Take a photo of the box markings.**
A photo of the markings (text) on the side of the box is required in case a part is needed for replacement. This helps our staff identify your product number to ensure you receive the correct parts.
- ⚠ **Take a photo of the damaged part (if applicable).**
A photo of the damage is always required to file a claim and get your replacement or refund processed quickly. Please make sure you have the box even if it is damaged.
- ⚠ **Send us an email with the images requested.**
Email us directly from marketplace where your item was purchased with the attached images and a description of your claim.

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