Foot Massager with Heat
Fußmassagegerät mit Wärme
Masseur de Pieds avec Chauffage
Masajeador de Pies con Calor
Massaggiatore per Piedi con Calore
Masażer do stóp z funkcją podgrzewania
Voetmassageapparaat met warmte

JS10033

USA office: Fontana **AUS** office: Truganina **ITA** office: Milano **POL** office: Gdańsk
DEU office: FDS GmbH, Neuer Höltingbaum 36, 22143 Hamburg, Deutschland **FRA** office: Saint Vigor d'Ymonville
GBR office: FDS Corporation Limited, Unit 4, Blackacre Road, Great Blakenham, Ipswich, Ip6 0FL, United Kingdom

If you're having difficulty, our friendly
customer team is always here to help.

USA:cs.us@costway.com
AUS:cs.au@costway.com
GBR:cs.uk@costway.com
ITA:cs.it@costway.com

POL:sklep@costway.com
DEU:cs.de@costway.com
FRA:cs.fr@costway.com
ESP:cs.es@costway.com

THIS INSTRUCTION BOOKLET CONTAINS **IMPORTANT** SAFETY INFORMATION.
PLEASE READ AND KEEP FOR FUTURE REFERENCE.

EN

Before You Start

- ⚠ Please read all instructions carefully.
- ⚠ Retain instructions for future reference.
- ⚠ Separate and count all parts and hardware.
- ⚠ Read through each step carefully and follow the proper order.
- ⚠ We recommend that, where possible, all items are assembled near to the area in which they will be placed in use, to avoid moving the product unnecessarily once assembled.
- ⚠ Always place the product on a flat, steady and stable surface.
- ⚠ Keep all small parts and packaging materials for this product away from babies and children as they potentially pose a serious choking hazard.

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EN


Technical Specifications:

Model	JS10033
Power Supply	24V, 50-60Hz
Rated Power	50W
Dimension	57 x 44.5 x 23 cm
Weight	8 KG
Package Include	1 x Massager
	1 x Power Adapter
	1 x Remote
	1 x User Manual

Massager Features:

- This electric foot massager provides a deep kneading shiatsu massage for feet, legs, ankles and calves.
- Heat therapy helps to accelerate blood circulation, maintain muscle and skin vitality.
- 5 Massage Modes: Adjustable air compression and massage styling.
- Remote control allows the elderly to enjoy massage easily.
- 360° adjustable stand.
- Electric foot massager integrates modern technology and ancient foot therapies to provide substantial health benefits besides stress relief. It improves the quality of life and strengthens the body.

Structure Diagram:

- A. Control Panel and Indicator
- B. Detachable and Washable Foot Cover
- C. Adjustable Stand
- D. Release Button
- E. Power Adapter
- F. Remote Control
- H. Vibration Motor

Cautions and Warnings General Safety:

Only use your foot massager as described in this manual.

Always unplug foot massager after using or before servicing (like detaching/washing the foot cover).

Do not operate your massager under blanket or pillow. Excessive heating may occur and cause fire, electric shock, or injury to you.

Do not place the device close to electric heaters, stoves or on place of direct sunlight.

Do not use the unit with wet feet. Keep your foot massager away from water, and wet or damp areas. Never place in water or liquid.

Do not use the device when the power cord or plug is damaged or soaked in water. Do not use device in moist or excessive dust environment to avoid damage.


This product is not recommended for children without parent's accompany.

For people whose feet are sensitive to the massaging nodes, this massager may be a little bit rough in the beginning. We suggest to lift your feet a bit and gradually your feet will adapt to the massager after several times use.


If you are heat-sensitive or experience discomfort while using the device, please turn off the unit.

This foot massager is not recommended for people with physical disabilities, sensory or mental disorders or lack of relevant experience.

Operation Guide:

1. On/Off: power button. Plug in the power cord, the massager starts in standby. The power indicator will flash.
2. Heat: heat therapy switch button.
3. Vibration: turn on or off the vibration function, and adjust 2 vibration options.
4. Intensity: Air compression adjust button. There are 5 different compression options from low to high.

5. Mode: 5 different mode options

Mode	Rolling/ Shiastu	Air Compression Intensity	Introduction & Suggestion	Heat and Vibration
Mode 1 	Forward/ reversal (alternative)	Default Gear 1 (low), adjustable	Highly recommend! Your feet will be moved horizontally by the air bags in the massager. The whole sole of the feet will be massaged. If you can wiggle knees lightly to adjust the feet angle to be in concert with the massager when it is working, you will enjoy more joyful massage.	Turn on or off based on personal preferences

Mode 2 	Forward only	Default Gear 3 (medium), adjustable	The rolling and shiastu will work forward only. There are more pauses in each circulation than Mode 1.	Turn on or off based on personal preferences
Mode 3 	Forward/ reversal (alternative)	Default Gear 3 (medium), adjustable	Your feet will not be moved horizontally by the air bags in the massager. The massage circulation will be shorter than Mode 1.	Turn on or off based on personal preferences
Mode 4 	None	Default Gear 3 (medium), adjustable	Only the air bags work to kneading your feet.	Turn on or off based on personal preferences
Mode 5 	Forward only	None, non-adjustable	No air compression intensity	Turn on or off based on personal preferences

Note: the massager will power off after 15 mins in usage.

Product Maintenance

- Make sure the machine is off and pull out the plug before cleaning.
- Please use a dry cloth to clean the surface of the product. Never use corrosive cleaning agents such as benzene or thinner.
- For the reusable foot cover, please make it dry before using it on the device.
- Please place the machine in a box and store in a safe, cool, dry and ventilated environment.
- Avoid scratching surface of the product by sharp objects.

Troubleshooting

Issue	Solution
The massager does not work	Ensure the unit is properly plugged into an appropriate outlet and turned on. If the power indicator lights up, but the massager does not work, please reach out to customer support team via order central.
The massager suddenly powers off	1. Ensure the unit is plugged in & powered on. 2. Time limit has been reached. The unit will be automatically power off after 15 minutes of usage. 3. The unit has overheated and the temperature protection feature was activated to shut off the unit. Wait until the unit has cooled down completely before resuming use. 4. Too much force has been applied to the massaging parts inside the unit. Make sure to keep your feet relaxed when inserted into the massager. Do not stand up while using the unit or exert undue force to be applied to the internal mechanisms of the unit.

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EN

Return / Damage Claim Instructions

- ⚠ **DO NOT discard the box / original packaging.**
In case a return is required, the item must be returned in original box. Without this your return will not be accepted.
- ⚠ **Take a photo of the box markings.**
A photo of the markings (text) on the side of the box is required in case a part is needed for replacement. This helps our staff identify your product number to ensure you receive the correct parts.
- ⚠ **Take a photo of the damaged part (if applicable).**
A photo of the damage is always required to file a claim and get your replacement or refund processed quickly. Please make sure you have the box even if it is damaged.
- ⚠ **Send us an email with the images requested.**
Email us directly from marketplace where your item was purchased with the attached images and a description of your claim.

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