



## TENS Unit Comparison Chart

Model Number	Pocket Pain Pro® PM400	Max Power Relief® PM500	Total Power + Heat™ PM800	Avail™ PM601
<b>Pain Modes</b>	3 Modes (Arm/shoulder,	6 Modes (Arm, Lower Back,	9 Modes (Arm, Lower Back,	6 Modes (Lower Back, Shoulder,
<b>Massage-like Modes</b>	2 Modes (Knead, Steady)	3 Modes (Tap, Knead, Rub)	3 Modes (Tap, Knead, Rub)	4 Modes (Steady, Knead, Tap,
<b>Intensity Levels</b>	10	15 (2x the power of the Pocket Pain Pro)	20	20
<b>Channels</b>	Single Channel	Single Channel	Single Channel	Dual Channel
<b>Bluetooth® App Controlled</b>				✓
<b>Wireless</b>				✓
<b>Heat</b>			✓	
<b>Pads</b>	1 Pair Long Life Pads™ (PMLLPAD)	1 Pair Long Life Pads™ (PMLLPAD)	2 Pair Gel Pads (PMGEL)	2 Wireless Pads (1 PMWPAD - M, 1
<b>Screen</b>		Easy-to-Read LCD	Easy-to-Read LCD	Smart phone
<b>Maximum Therapy Duration</b>	15 Min	15 Min	30 Min	60 Min
<b>Warranty</b>	1 year limited warranty – unit 30 day limited warranty – Long Life Pads™	1 year limited warranty – unit 30 day limited warranty – Long Life Pads	1 year limited warranty – unit 30 day limited warranty – Gel Refills	1 year limited warranty – unit 30 day limited warranty – Wireless Pads
<b>Battery</b>	2 AAA batteries included	2 AAA batteries included	Rechargeable Lithium Ion Battery	Rechargeable Lithium Ion Battery
<b>Belt Clip</b>		✓		
<b>QSG &amp; Pad Placement Guide</b>	✓	✓	✓	Within App
<b>Pad Holder</b>	✓	✓	✓	✓
<b>Travel Case</b>		✓		✓