

# CERAGEM

MY FUTURE HOME



## CES 2025 Six Innovation Awards

### Health & Wellness Products

#### CES 2025 INNOVATION AWARD PRODUCT

- **DIGITAL HEALTH (4)** - Home Therapy Booth / Home Medical Bed / Smart Headboard / AI Medi Water
- **FITNESS (1)** - Therapeutic Spa Bed
- **HUMAN SECURITY FOR ALL (1)** - Neuro Wellness Enhancer



# CERAGEM V7

## COLOR VARIANTS



Black

White


Brown

**Consumption Power:** 480 VA

**Working Voltage:** AC 220V, 50/60 Hz

**Unit Weight:** 162 lbs **w/ Box:** 166 lbs

**Tolerance Weight:** 300lbs



## DIMENSIONS

### BED

width	28.5"
height	17"
depth (bed opened)	80.5"
depth (bed folded)	47.7"

### BOX

width	32"
height	19"
depth	53"

## INCLUDED ATTACHMENTS



Abdominal Thermal Projector



UPGRADED  
Air-Cell lower  
body massager



Class II Cleared

## KEY POINTS

PATENT

- Auto Spine Scanning Technology
  - Automatically scans user's spinal curvature, length & weight distribution

- Mode 8 (cervical)

PATENT

- Tilting massage rollers for deeper neck massage

- Tilting massage rollers for deeper neck massage

Voice guidance is provided to support effective cervical massage.

- The internal applicator is positioned at the T10 (10th thoracic vertebra) to encourage gentle head rotation, helping stretch the back, chest muscles, and shoulders.

## PRICE

**MSRP:** \$7,999 **White Glove:** \$200

## RETURN POLICY

- return within 14 days of delivery date
- 20% restocking fee

## WARRANTY POLICY

- Standard Warranty** (Year 1): Parts & Labor
- Extended Coverage** (Years 2-3): Parts Only
- Additional Support:** 10-year manufacturer parts availability (main PCB & motor)

\*Covered components: Rollers, heating systems, frame structures, display units & electronics

## V7 REDEFINES IN-HOME WELLNESS

With minimalist design & intuitive function, the V7 fits seamlessly into the heart of your home- whether it's a cozy bedroom, a calm living space, or even a child's room for gentle, supervised use.



## WELLNESS THAT LIVES WHERE YOU DO- DESIGNED FOR COMFORT, CONNECTION & SHARED CARE

- Elegant enough for a living room
- Quiet & sleek for the bedroom
- Gentle & safe for family-friendly environments
- Built-in programs suited for kids, parents & grandparents
- Refined curves & a calming color palette that promotes relaxation



## MODERN WELLNESS MEETS TIMELESS DESIGN

The V7 blends effortlessly into any living space with its slim profile, neutral tones, and understated elegance. No longer hidden away, this is a wellness experience you'll be proud to display.



# CERAGEM V7 WHATS NEW



## Mode 8 Cervical (neck/shoulders)

- Allows users to rotate their heads in response to voice guidance actively
- The head rotation feature during the massage, the V7 can now reach side neck muscles, like the obliquus capitis superior and sternocleidomastoid, that the V6 couldn't access before.
- Provides a more intense & effective pressure to targeted areas
- **Speed adjustment not available; temperature cannot exceed 113°F**



## Simplified Massage Modes

- Focus Mode → Full Spine Mode
- Manual Mode → Focused Heat Mode
- Semi-Auto Mode → Focused Massage Mode
- Intuitive mode names



## DIRECT HEATING MASSAGE BALL



### Point 1

- **Fast heat-up speed**  
Quickly reaches up to 65°C in just 10 minutes.

### Point 2

- **Minimal heat loss**  
By filling the gap between the inner heated part and the insert part with the rotating pin, heat loss is minimized

### Point 3

- **Soft touch with silicone cushion**  
Finished with silicone material, providing a soft massage feel


## 4<sup>th</sup> GENERATION

- V6 & V7
- **Dual-row jade heads**
  - closer spinal contact & curvature-tracing, providing more precise, personalized massage
- Scans the user's spine
- Measures the curvature
- **Analyzes weight distribution**
  - Analyzes the user's body weight and ensures that all intensity levels always feel the same
- Spinal thermal ceramic heads massage vertically, horizontally & **tilt**



## Upgraded Aircell Leg Massager Intensity

- More precise & gentle compression
- More detailed control
- 3 massage modes
  - feet to calves
  - whole legs
  - calves to thighs
- 18 mins



## Tilting Massage Rollers

- Massage rollers tilt to match spinal curvature for close contact massage.
- Closer & precise contact around the neck, shoulder & upper spine
- Split rollers (front & rear) give highly curved areas like neck & lower back, w/o lifting or gaps.




## Comfort Mat & Soft Touch




- Thicker foam & dual inner fabric layers for a more comfortable massage
- Softer pelvic massage experience

# CERAGEM V7

## WHATS INCLUDED

## ATTACHMENTS



 <p><b>Aircell Leg Massager</b></p>	<p>Helps blood circulation using segmented strength control (limbs compression &amp; circulation device)</p> <ul style="list-style-type: none"><li>• Improves circulation &amp; alleviates tense legs</li><li>• Repeatedly expands &amp; contracts</li></ul>
 <p><b>Abdominal Projector</b></p>	<ul style="list-style-type: none"><li>• Heats up to 140F°</li><li>• soothes abdomen</li></ul>

 <p><b>Main Body</b></p>	 <p><b>Abdominal Projector</b></p>
 <p><b>Remote Control</b></p>	 <p><b>Outer Cloth</b></p>
 <p><b>Head Cushion</b></p>	 <p><b>Aircell Leg Massager</b></p>
 <p><b>User Manual</b></p>	 <p><b>Power Cord</b></p>
 <p><b>Installation Guide</b></p>	 <p><b>Abdominal Projector Cover</b></p>

## V7 CHARACTERISTICS

 <p><b>3D Motion Spinal Thermal Projector</b></p>	<p>Reinforcement of pressure for the neck + shoulders &amp; close massage along the spine (pressure stimulation)</p>
 <p><b>Customized Massage</b></p>	<p>Functions tailored to the user's body type. Customized massage tailored to each user's body type by measuring the length &amp; curvature of the spine.</p>
 <p><b>3 Step Speed Control</b></p>	<p>Adjusting the massage speed in 3 steps, as fast or slow as the user desires or requires.</p>
 <p><b>Auto Start</b></p>	<p>Using the weight sensor, massage begins automatically after the user lies down while the power is on</p>
 <p><b>Easy Sliding</b></p>	<p>Mat structure is easily opened &amp; folded with wheels incorporating the actuator technology</p>

# V7 CONTROLLER + CONNECTOR



Find & select Master V7-XX from bluetooth device search list.


When the connection is made, a bluetooth icon will appear at the top of the screen

**Please note:** USB port has been removed

# V7 MESSAGE MODES

## STANDARD MODES

Mode 1 (Standard)	heat + massage neck, back, waist & pelvis (muscle pain relief)	36m 30s
Mode 2 (Relax)	light & gentle massage back, waist, and pelvis	18 mins
Mode 3 (Special)	#1 w/ intensified heat + massage (pain relief)	75 mins
Mode 4 (Stable)	focused massage + heat on back gentle massage on waist, pelvis & neck	36m 30s
Mode 5 (Circulation)	similar to #4 w/ slightly different massage pattern	36m 30s
Mode 6 (Study)	Focused heat + massage neck, back, and pelvis (ideal after prolonged sitting)	36m 30s
Mode 7 (Pressure Therapy)	assists w/ herniated discs & spinal stenosis; traction intensity automatically adjusts	36m 30s



## DAILY ROUTINE MODES

Mode 11 (Morning/ Stretch)	Loosens stiff neck, back & waist from overnight tension	36m 30s
Mode 12 (Evening)	Releases daily neck & back muscle tension	36m 30s

## DAILY LIFESTYLE MODES

Mode 13 (Office Worker/ Driver)	focuses on relieving back & waist fatigue from long work or driving hours (V6-> Energy Mode)	36m 30s
Mode 14 (Housework / Childcare)	eases back & waist tension after chores & childcare	36m 30s
Mode 15 (Sports/Golf)	relieves back, pelvis & waist tension after physical activities	36m 30s

## SELECTABLE MODES

Full Spine	massage for the entire spine or a specific vertebra	90 mins
Focused Heating	provides heat to selected spinal area	60 mins
Focused Heating (Master)	advanced focused heating mode with more control	60 mins
Focused Massage	massage for a selected spinal area	60 mins
Focused Massage (Master)	advanced version of focused massage mode.	60 mins



## INTENSIVE MODE FOR EACH PART




Mode 8 (Neck/Shoulders)	focuses on neck + upper back; interactive; limited temperature & intensity	36m 30s
Mode 9 (Upper back)	heat + massage on upper back; pelvis & neck massage	18 mins
Mode 10 (Lower Back)	intense massage + heat on lower back, waist & pelvis	36m 30s

new semi / manual modes


# V7 SCREEN DISPLAY



## TEMPERATURE SETUP



## AIRCELL SETUP



Wait!

- When the supporting mat is closed, an image of the closed supporting mat will be displayed as the product image.
- The screen of the remote control screen is subject to change without notice for performance improvement.

# SAFETY PRECAUTIONS



## WARNING

To avoid risk or injury:



- Those with sensitive skin must use the product at lower temperatures & be cautious when using for an extended period, even if it is set at a low temperature, due to the risk of burning
- Do not use for more than 2-3 mins on any one part of the body when setting a temperature over 113°F
- Do not use the product if you have recently undergone skin-related surgery
- Do not use the product if you have recently undergone surgery for varicose veins, etc.
- Stop using the product if the pain gets worse during use
- Do not use if there is bleeding in the area to which you intend to apply it, or if you are experiencing pain or inflammation of an unknown cause
- If the pressure of the air cell massager differs significantly, call customer service
- Do not use the air cell on bare skin
- If using the product on the abdomen, do not use it on the ribs
- Do not use on bare skin
- Do not use if you have had any major surgeries, including back & hip
- Stop using the product if you start feeling sick or uncomfortable, or if the massage seems too painful
- Do not sleep while using the product
- Do not use the product after drinking alcohol or sleep
- Do not let kids or pets play around the product during operation



## Consult Doctor



- Pregnant women
- Osteoporosis
- Poor physical health / undergoing any medical treatment
- Higher body temperature. IE: Chills
- Symptoms of an acute skin disorder, malignant tumor, or other skin disorders
- Suffering from back, neck, shoulder, or hip pain
- Have a pacemaker, defibrillator, other artificial implants, or other medical devices
- Sensitive reactions or allergy to heat
- Deep vein thrombosis, blood clot, or edema
- Inflammation of the skin
- Thrombophlebitis
- Ischemic vascular disease
- Severe peripheral neuropathy
- Have any external wounds or infections on the skin where the product will come in direct contact
- No evidence suggests the product interferes w/ or influences medications, but they should still consult their doctor
- Suffer from pain of an unidentified cause
- Suffer from weak circulation in legs or heart disease



## Must Be Under Supervision



- reduced physical, sensory, or mental capabilities
- Lacks appropriate experience or training
- children

# CERAGEM V5

## TRAINING GUIDE

### COLOR VARIANTS



Brown



White



Black

**Consumption Power:** 380 VA  
**Working Voltage:** 100-127V-, 50/60 Hz  
**Unit Weight:** 109lbs (± 6lbs)  
**Weight w/ Box:** 136lbs (± 6lbs)  
**Tolerance Weight:** 298lbs



### DIMENSIONS

#### BED (folded)

width	27.6"
height	17.7"
depth	50.4"

#### BED (open)

width	27.6"
height	17.7"
depth	80.3"

#### BOX

width	53"
height	20"
depth	31"

### INCLUDED ATTACHMENTS



Abdominal Vibration Projector

### PRICE

**MSRP:** \$4,599 **White Glove:** \$200

### RETURN POLICY

- return within 14 days of delivery date
- 20% restocking fee

### WARRANTY POLICY

- **Year 1:** Parts & Labor
- **Year 2:** Parts Only
- **Remaining 7 years:** Only CeraCore is covered
- **Labor not included beyond year 1**

### KEY POINTS

- PATENT** • Auto Spinal Scan
- Thermal & vibrating abdominal projector
- PATENT** • 3rd generation CeraCore engine
- PATENT** • Sliding leg design
- 17 pre-programmed massage modes
  - 12 automatic + 5 Manual

# MASTER V5 Health Benefits

## Thermal & Infrared Heat Therapy

- Relieves muscle tension & stiffness
- Improves local blood circulation where heat is applied
- Helps with fatigue recovery by relaxing deep muscle layers
- May support immune function through improved circulation & detoxification

## Spinal Stretching & Acupressure

- Promotes spinal alignment & posture correction
- Helps decompress vertebrae to reduce back & neck pressure
- Targets key acupressure points for enhanced relaxation & body balance
- May support relief of sciatica, disc-related discomfort, or spinal fatigue



# MASTER V5 Health Benefits

## Stress & Mental Wellness Support

- Encourages mental relaxation through Ceragem Sound Therapy (nature, classical, functional music)
- Supports meditative breathing and mindfulness through dedicated modes
- Reduces stress-related muscle tension, especially in the neck & shoulders

## FDA-Cleared Benefits

- Temporary relief of minor muscle and joint pain and stiffness
- Temporary relief of minor joint pain associated with arthritis
- Temporary increase in local circulation
- Relaxation of muscles

## Core & Abdominal Stimulation

- Helps stimulate deep core muscles
- Supports digestive comfort & abdominal relaxation
- Promotes a sense of full-body balance & comfort



# CERAGEM V5



## WHAT MAKES V5 FDA II CLEARED?

Provides muscle relaxation therapy by delivering heat and soothing massage.

### Provides topical heat for:

- Temporary relief of minor muscle & joint pain stiffness
- Temporary relief of minor joint pain associated with arthritis
- Temporary increase in local circulation where applied
- Relaxation of muscles

## V5 CHARACTERISTICS

 <b>Body Sensing System</b>	Automatically measures spinal length, curvature, and weight distribution to personalize the massage path.
<b>17</b> <b>17 Core Massage Modes</b>	Includes automatic, manual, semi-auto, and targeted master modes (C1-L5) for personalized care.
 <b>Master Mode</b>	Allows users to select specific vertebrae (C1-L5) for focused thermal massage based on their needs.
 <b>Built-In Bluetooth Speaker</b>	Pairs with smart devices for personal audio
 <b>Abdominal Vibration Projector</b>	Provides vibration and heat in <b>3 modes</b> with <b>2 intensity</b> levels – supports core relaxation and comfort.
 <b>Patented Sliding Design</b>	Leg portion retracts for space-saving; doubles as a heated ottoman when not in use
<b>3rd gen</b> <b>3rd Generation CeraCore Engine</b>	Delivers vertical and horizontal spinal traction
 <b>Health-terior</b>	Designed for interior harmony- Perfect for bedrooms, living rooms & small spaces
 <b>Voice-Guided Interface</b>	Step-by-step voice instructions assist with operation and program selection

# CERAGEM V5

## THERMAL MASSAGE ROLLERS

Massage rollers are made of **Ceratonic**, a special material created by Ceragem, made with **jade, yellow clay, maekbanseok** (mineral stone used in traditional Korean wellness), etc. These materials are beneficial to the body.


## Ceracore Engine 3<sup>rd</sup> Generation




### 3<sup>RD</sup> GENERATION

- V3, V4 & V5
- Scans user's spine
- Measures the curvature
- Spinal thermal ceramic heads massage vertically & horizontally

# V5 CONTROLLER + CONNECTOR




# V5 REMOTE CONTROL




- Main screen -


1	product no.
2	Sound source (Ceragem sound / USB / Bluetooth / Voice guidance)
3	Sound pause
4	Volume (0~20)
5	Time remaining (in minutes)
6	Projector location
7	Mode (Automatic mode A,1~11/Intensive/Semi-automatic/Manual)
8	Strength display (1-9)
9	Internal projector temperature (86~149°F)
10	Display of connection of the abdominal vibration projector



When using the thermal sofa



When setting the temperature



When setting the vibration projector

Wait!

- When the supporting mat is closed, the "thermal sofa" icon will be displayed.
- When 10 seconds have passed after the remote control was last operated, the screen will become dark. (When a button is operated, the screen will become bright again.)

# V5 MESSAGE MODES

## STANDARD MODES

Mode & Program Name	Description	Focus Area	Duration
Mode A Standard	Full spine care mode with thermal massage across 13 key spinal points	Full Body	30m 30s
Mode 1 Special	Combines full spinal care with Sleep Mode	Full Body	75m
Mode 2 Relax	Provides fast, comfortable muscle relaxation	Back, Waist, Pelvis	18m
Mode 3 Stable	Focused care on back area, especially around the navel	Back Only	36m 30s
Mode 4 Circulation	Focused massage on the upper back/ chest area	Back Only	36m 30s
Mode 5 Basic Health	Targets concentrated thermal stimulation on thoracic and lumbar	Back Only	36m 30s
Mode 6 Energy	Focused massage & thermal stimulation on the lumbar region	Lower Back, Pelvis	36m 30s
Mode 7 Study	Targeted care for cervical spine	Neck Area	36m 30s

Automatic mode	
⊖ A	ModeA(Standard)
⊖ 1	Mode1(Special)
⊖ 2	Mode2(Relax)
⊖ 3	Mode3(Stable)
⊖ 4	Mode4(Circulation)
⊖ 5	Mode5(Basic health)
⊖ 6	Mode6(Energy)
⊖ 7	Mode7(Study)
⊖ 8	Mode8(Rest)
⊖ 9	Mode9(Body line)
⊖ 10	Mode10 (Meditation-relax)
⊖ 11	Mode11 (Meditation-stable)
Selection mode	
👤	Manual mode
👤👤	Intensive mode
👤👤👤	Semi-automatic mode

Back label

## STANDARD MODES

Mode & Program Name	Description	Focus Area	Duration
Mode 8 Rest (Sleep)	Recommended for those seeking comfortable sleep	Neck, Back	36m 30s
Mode 9 Body Line	Feels like a fully-body massage, enhancing body alignment	Back, Waist, Pelvis	36m 30s
Mode 10 Meditation - Relax	A short program paired with body scan meditation for light relaxation	Back, Waist, Pelvis	18m
Mode 11 Meditation - Stable	Promotes mental & physical calm through body scan meditation	Back Only	36m 30s

## SELECTABLE MODES

Manual Mode	Move the rollers freely to your desired position (with heat)	60m
Manual Mode (Master)	Input vertebra number & receive targeted manual massage at that location	60m
Intensive Mode (Focused)	massage for the entire spine or a specific vertebra	90m
Semi-Automatic	Move the rollers to a chosen position (repeats up/down + heat)	60m
Semi-Automatic (Master)	Input vertebra number & receive targeted massage	60m

# SAFETY PRECAUTIONS

## **Who Should Not Use The Device Without Medical Advice**

- Pregnant individuals
- Those with malignancies, high fever, or severe osteoporosis
- People recovering from back, neck, shoulder, or hip surgery
- Anyone with a pacemaker, defibrillator, or medical implant
- Users with heat sensitivity or skin conditions
- Children or individuals unable to adjust temperature independently

## **General Warnings**

- Do not sleep while using the product
- Do not use the device after consuming alcohol
- Keep children & pets away during operation
- Never disassemble, repair, or modify the product yourself
- Avoid using in wet, humid, or high-dust environments

## **Heat-Related Cautions**

- Always wear long-sleeved clothing during use
- Do not use the abdominal vibration projector or external projector on the same spot for more than 2-3 minutes.
- Burn risk: Be cautious with high heat settings or prolonged contact with sensitive areas.
- Never place external projectors on the eyes or face
- Use the external projector only with the Ceragem-provided cover.

## **Electrical Safety**

- Unplug the device if not in use or during electrical storms
- Do not operate with a damaged cord or plug
- Never pull the cord to unplug or touch the plug with wet hands
- Ensure the power source matches the 100-127V~ rated voltage

## **Usage Guidelines**

- Use 1-2 times per day as recommended
- Do not use while lying face down, except for designated abdominal modes
- Keep the remote & cables secure- don't yank or throw them
- Sit and lie down carefully to avoid falls or injury during use