



TENS Accessories Comparison Chart

Model Number	Long Life Pads™ PMLLPAD	Long Life Pads™ - Large PMLLPAD - L	Gel Refills PMGEL	Avail Wireless Pad - Medium PMWPAD - M	Avail Wireless Pad - Large PMWPAD - L
Reusable	Lab tested for up to 150 uses	Lab tested for up to 150 uses	Lab tested for up to 30 uses	Lab tested for up to 30 uses	Lab tested for up to 30 uses
Washable	Yes, up to 10 times (use tap water, air dry)	Yes, up to 10 times (use tap water, air dry)	No	Yes, up to 2 times (use moistened cloth, do not use paper cloth)	Yes, up to 2 times (use moistened cloth, do not use paper cloth)
Size	4" x 2.5" each 1 Pair	4.5" x 3" each 1 Pair	1.6" x 1.6" each 3 Pair	3.1" x 7" each 1 Pad	3.3" x 8.6" each 1 Pad
Pre-gelled, Self-adhesive	✓	✓	✓	✓	✓
Soft, flexible material	✓	✓	✓	✓	✓
Single Person Use	Yes	Yes	Yes	Yes	Yes
Compatible With	PM3029, PM3030, PM3031, PM3032, PM400	PM3029, PM3030, PM3031, PM3032, PM400	PM311	PM601	PM601
Pad/Gel Subscription Discount	✓	✓	✓	✓	✓
Pad Holder	✓	✓			
Warranty	30 day limited warranty	30 day limited warranty	30 day limited warranty	30 day limited warranty	30 day limited warranty

5/8/2019